

Monoclonal Antibody to Human Chorionic Gonadotropin hCG (A2010003-4D9)

| | |
|-----------------------------|---|
| Catalog number: | A2010003-4D9 |
| Lot number: | 0119A0003-4D9 |
| Expiration date: | 03/10/2023 |
| Others Name(s): | Mab to hCG |
| Host: | Sheep |
| Clone: | 4D9 |
| Isotype: | IgG1 |
| Purity: | >90% (protein A column) |
| Concentration: | 1.0 mg/mL |
| Storage: | -20°C. Avoid repeated freeze/thaw cycles. |
| Supplied as: | Frozen solution |
| Buffer: | PBS 1x + 0.05% sodium azide |
| Keywords: | hCG antibody, Human Chorionic Gonadotropin |
| Product Description: | Monoclonal antibody to human chorionic gonadotropin hCG. Antibody was obtained by immunizing sheep with human hCG (whole protein). Epitope mapping shows that the antibody recognizes both intact hCG and hCG beta subunit. Two clones are available (4D9 and 7H2). |

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.