

Monoclonal Antibody to Phenytoin (A2010005-1E5)

Catalog number:	A2010005-1E5
Lot number:	0119A0005-1E5
Expiration date:	03/13/2023
Others Name(s):	Mab to Phenytoin
Host:	Mouse
Clone:	1E5
Isotype:	IgG1
Purity:	>90% (protein A column)
Concentration:	1.0 mg/mL
Storage:	-20°C. Avoid repeated freeze/thaw cycles.
Supplied as:	Frozen solution
Buffer:	PBS 1x + 0.05% sodium azide
Keywords:	Phenytoin, antibody
Product Description:	Mouse monoclonal antibody to phenytoin (a drug used to treat seizures, heart arrhythmia, and neuropathic pain). Six clones are available: 1E5, 4D8, 8H9, 2C6, 5H7 and 9F3. Antibodies can be used in ELISA formats as a well as western blot and other immunodetection methods.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.