

Human Myoglobin Antibody

Catalog number:	A2010008
Lot number:	0119A0008
Expiration date:	03/27/2023
Package Size:	1 mg
Purity:	Affinity purified
Host:	Goat
Antigen:	Purified human myoglobin from human heart
Clonality:	Polyclonal
Concentration:	4 mg/mL
Buffer:	Tris pH 7.4
Preservative:	0.05% sodium azide
Storage:	2-8°C (short term). -20°C (long term). Avoid repeated freeze-thaw.
Product Description:	Myoglobin is an oxygen-binding protein found in muscle cells, including heart muscle cells. This polyclonal antibody targets specifically human myoglobin found in heart cells. Presence of myoglobin in serum is often indicative to injury to the heart muscle.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.