

Digoxigenin antibody

Catalog number:	A2010019
Lot number:	0119A0019
Expiration date:	04/06/2023
Package Size:	0.25 mg
Concentration:	1 mg/mL
Other names:	antibody to Digoxigenin, anti-Digoxigenin antibody
Supplied as:	Solution (PBS1X)
Host:	Goat
Clonality:	Polyclonal
Isotype:	IgG
Purification:	Affinity purified
Preservative:	0.05% sodium azide
Storage:	2-8°C (short term), -20°C (long term). Avoid repeated freeze-thaw.
Product Description:	Digoxigenin is a steroid molecule extracted from the flowers and leaves of <i>Digitalis purpurea</i> , <i>Digitalis orientalis</i> and <i>Digitalis lanata</i> . Digoxigenin reacts with sugars to form glycosides. This polyclonal antibody to Digoxigenin can be used in ELISA and in situ hybridization experiments.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.