

## Monoclonal Antibody to Digoxin

<b>Catalog number:</b>	A2010020
<b>Lot number:</b>	0119A0020
<b>Expiration date:</b>	04/06/2023
<b>Package Size:</b>	0.2 mg
<b>Other names:</b>	Antibody to digoxin
<b>Supplied as:</b>	Liquid Solution
<b>Concentration:</b>	2 mg/mL
<b>Host:</b>	Mouse
<b>Antigen:</b>	Digoxin molecule conjugated to HSA
<b>Clonality:</b>	Monoclonal (Clone 3H9)
<b>Preservative:</b>	0.05% sodium azide
<b>Storage:</b>	2-8°C (short term), -20°C (long term). Avoid repeated freeze-thaw.
<b>Specificity:</b>	Lanatoside C 99% beta-methyl digoxin 98% Deslanoside 97% Digoxigenin 44% Digotoxin 15% Quabain 1.6%
<b>Product Description:</b>	Monoclonal antibody to Digoxin, a drug used as a beta-blocker in the treatment of several medical conditions such as congestive heart failure (CHF).

---

**WARNING:** This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

**SAFETY:** User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

**WARRANTY AND LIMITATIONS:** Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.