

Carbamazepine Antibody

Catalog number:	A2010024
Lot number:	0119A0024
Expiration date:	4/15/23
Package Size:	0.25 mg
Other names:	Antibody to Carbamazepine
Supplied as:	Liquid Solution
Concentration:	4 mg/mL
Type:	Ig
Host:	Sheep
Specificity:	Carbamazepine (100%), Carbamazepine 10,11-epoxide (11%)
Cross-reaction:	Carbamazepine 10,11-epoxide (11%), 10,11-Dihydro-10-hydroxycarbamazepine (1%)
Purification:	Purified
Immunogen	Carbamazepine-BSA
Preservative:	0.05% sodium azide
Storage:	2-8°C (short term), -20°C (long term). Avoid repeated freeze-thaw.

Product Description: Carbamazepine is an anti-seizure drug belonging to the class of anti-epileptics and anti-convulsants. This sheep polyclonal antibody is specific to carbamazepine and does not significantly cross-react with carbamazepine-derived metabolites.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.