

L-Thyroxine

Catalog number:	C2010004
Lot number:	0119C0004
Expiration date:	04/06/2023
Molecular Weight:	776.9 g/mol
Purity:	> 98%
Package Size:	100 mg
Appearance:	White to Light Yellow Powder
Storage:	-20°. Avoid repeated freeze-thaw.
Product Description:	L-Thyroxine (T ₄) is an iodine-containing hormone produced by the thyroid glands. T ₄ plays a central role in cellular physiology as it regulates both growth and development by acting directly and indirectly on the level of expression of various enzymes and proteins.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.