

Recombinant Human Serum Amyloid A1 (SAA)

Catalog number:	P2010014
Lot number:	0119P0014
Expiration date:	4/21/23
Package Size:	1.0 mg and bulk
Other names:	SAA human
Supplied as:	Lyophilized powder
Molecular Weight:	12 kDa
Purity:	> 95% pure
Buffer:	Reconstitute in PBS at 0.1 mg/mL.
Storage:	2-8°C (short term), -20°C (long term). Avoid repeated freeze-thaw.

Product Description:

Human Serum Amyloid A1 (SAA) produced in *E. coli*. SAA is an apolipoprotein produced in the liver during cytokine-based responses to inflammation. SAA levels in the blood can be used to monitor acute and chronic inflammation.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.