

Folate Binding Protein, Bovine

Catalog number:	P2010016
Lot number:	0119P0016
Expiration date:	4/21/23
Package Size:	0.5 mg
Other names:	FBP from Bovine Milk
Supplied as:	2 mg/mL solution
Molecular Weight:	26 kDa
Purity:	> 95% pure
Buffer:	PBS
Storage:	2-8°C (short term), -20°C (long term). Avoid repeated freeze-thaw.

Product Description:

Folate Binding Protein (FBP) is abundant protein bound in milk as a monomer. The protein contains 8 disulfides bonds and is glycosylated. Its main role is to bind folate. As such, FBP plays a crucial role in the metabolism of folate.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.