

SAA Antibody Polyclonal

Catalog number:	A2010028
Lot number:	0119A0028
Expiration date:	04/21/2023
Others Name(s):	Polyclonal Antibody to Human Serum Amyloid Protein A1
Host:	Rabbit
Clone:	Polyclonal
Isotype:	IgG
Purity:	>90% (protein A column)
Concentration:	5 mg/mL
Storage:	-20°C. Avoid repeated freeze/thaw cycles.
Supplied as:	Frozen solution
Buffer:	PBS 1x + 0.05% sodium azide
Keywords:	SAA, Anti-Serum amyloid A
Product Description:	Polyclonal antibody to human Anti-Serum amyloid A. This rabbit antibody was obtained by using human SAA as immunogen. This antibody has been validated for use in ELISA, WB and homogenous immunoassays. BSA cannot be used as a blocking agent with this antibody.

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes. Furthermore, as a universal safety precaution, we recommend that this product be treated as hazardous until information to the contrary becomes available.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.