


## 25-OH Vitamin D BSA Conjugate

<b>Catalog number:</b>	P2010018
<b>Lot number:</b>	0119P0018
<b>Expiration date:</b>	07/22/2023
<b>Others Name(s):</b>	BSA-Vitamin D, Vitamin D-BSA, BSA conjugated to 25-OH Vitamin D, Vitamin D Bovine Serum Albumin Conjugate, 25-OH Vitamin D [BSA].
<b>Molecular Weight:</b>	70 kDa (average)
<b>Concentration:</b>	1 mg/mL
<b>Amount:</b>	0.5 mg
<b>Purity:</b>	>95%
<b>Storage:</b>	-20°C. Avoid repeated freeze/thaw cycles.
<b>Supplied as:</b>	Liquid solution
<b>Buffer:</b>	1X Phosphate Buffer Saline, pH 7.4.
<b>Keywords:</b>	Vitamin D BSA Conjugate, CALCIDIOL; CALCIFEDIOL; CALCIFEDIOL MONOHYDRATE; 25-OH-D3; 25HYDROXYVITAMIN D3; 25-hydroxyvitamin d3 monohydrate; 25-HYDROXYCHOLECALCIFEROL; 25-hydroxycholecalciferol monohydrate; (5z, 7e)-9, 10-secocholesta-5, 7, 10(19)-triene-3b, 25-diol monohydrate; VITAMIN D3, 25-HYDROXY CALCIFEDIOL 1-HYDRATE; 25-OH-Vitamin D3

**Product Description:** Concentrated solution of Bovine Serum Albumin (BSA) conjugated to 25-OH Vitamin D3 (a pre-hormone produced in the liver by hydroxylation of Vitamin D3). A carbon spacer arm links BSA to the 25-OH Vitamin D3 moiety. Each BSA molecule contains several 25-OH Vitamin D3 units. Dilutions of the provided conjugate can be made in 1X Phosphate Buffer Saline, pH 7.4.



25-OH Vitamin D

**WARNING:** This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

**SAFETY:** User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes.

**WARRANTY AND LIMITATIONS:** Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.