


Everolimus BSA Conjugate

Catalog number:	P2010024
Lot number:	0120P0024
Expiration date:	2024-04
Others Name(s):	BSA-Everolimus, Everolimus-BSA, BSA conjugated to Everolimus, 40-O-(2-hydroxyethyl)-sirolimus BSA, Everolimus Bovine Serum Albumin Conjugate, Everolimus [BSA].
Molecular Weight:	71 kDa (average)
Concentration:	1 mg/mL
Amount:	0.5 mg
Purity:	>95%
Storage:	-20°C. Avoid repeated freeze/thaw cycles.
Supplied as:	Liquid solution
Buffer:	1X Phosphate Buffer Saline, pH 7.4.
Keywords:	BSA conjugate, Everolimus, Certican, Zortress, Afinitor, Votubia, Evertor

Product Description: Concentrated solution of Bovine Serum Albumin (BSA) conjugated to Everolimus. A carbon spacer arm of 19.0 Å links BSA to the Everolimus moiety. Each BSA molecule contains several Everolimus units. Dilutions of the provided conjugate can be made in 1X Phosphate Buffer Saline, pH 7.4.



Structure of Everolimus

WARNING: This product is for **Research Use Only (RUO)**. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.