



Intrinsic Factor Monoclonal Antibody

Catalog number:	A2010044
Lot number:	0120A0044
Expiration date:	2024-04
Clonality:	Monoclonal (Clone 17A9)
Immunogen:	Intrinsic Factor from porcine gastric mucosa
Host:	Mouse
Concentration:	5 mg/mL
Amount:	1 mg
Purity:	> 95% (FPLC)
Isotype:	IgG1
Storage:	2-8°C.
Supplied as:	Liquid solution
Buffer:	50 mM Tris-HCl pH 7.0 +500 mM NaCl + 0.05% sodium azide

Product Description:

Mouse monoclonal antibody to the Intrinsic Factor (IF) from porcine gastric mucosa. Highly pure and highly concentrated, this antibody has several applications that include immunoassay development, ELISA, WB and immune-histochemistry.

WARNING: This product is for **Research Use Only (RUO)** or Manufacturing. This product is not for administration to humans or animals. This product is not for human or veterinary diagnostic or therapeutic use.

SAFETY: User must review the **Material Safety Data Sheet (MSDS)** provided with this product. In addition, this product should not be ingested, inhaled, swallowed or put in contact with eyes, skin or clothes.

WARRANTY AND LIMITATIONS: Molecular Depot makes no warranty or guarantee of any kind (written, oral, expressed or implied) of the fitness or suitability of this product to any user's particular purpose or application. Furthermore, in no event will Molecular Depot LLC have any obligation or liability for any direct, indirect, incidental or consequential use of this product. All sales of this product are final and non-refundable.